

# PapooSwitch

Better than a classic external switch, PapooSwitch is an ultrasensitive tactile button



## Why choosing PapooSwitch?

### PapooSwitch. **Tactile.**

Sensitive to the approach of the human body, PapooSwitch acts like a switch when brushed.



### PapooSwitch. **Universal.**

It drives any electronic system that has an external switch input, like Papoo.



### PapooSwitch. **Adjustable.**

Sensitiveness is adjustable thanks to buttons located under your PapooSwitch. You can activate the bright signal and the tone to make the detection noticeable.



### Technical specifications



USB version also available !  
PapooSwitch USB then acts like the left or right click of your mouse.






Join Papoo's world and discover the whole range  
[www.papoo.fr](http://www.papoo.fr) / +33(0)466281146

# PapooSwitch



Reliable. Simple. Compatible. Adjustable to any situation.

## CASE n°1





Thanks to its hypersensitivity, you activate the PapooSwitch with any body part. You just need to approach PapooSwitch surface up to 1 centimeter.

## CASE n°2



PapooSwitch is adjustable on any surface and in any mobility condition.

## CASE n°3



PapooSwitch can also be used with a USB to be connected to your computer for instance. Therefore, it plays the role of the right or left click of a mouse.